






















































## Common screw/screw driver systems of AS

Brand	System	Level	Platform diameter	Type of screw head (Straight ONLY)		Pictures of Straight Screw Hand Driver	ANGLED SCREW DRIVER:	Original Custom angled screw driver ref code:	Picture of ANGLED SCREW DRIVER					
<b>3i Biomet</b>	Certain	Fixture	3.4		1.20mm		✓	ASCDT30 (30mm)						
			4.1											
			5.0											
	External Hex	Fixture	3.4											
			4.1											
			5.0											
<b>Dentsply</b>	AstraTech OsseoSpeed TX	Fixture	3.0		1.27mm		✓	P-DYNDRIVER						
			3.5/4.0											
			4.5/5.0											
	UniAbutment	Abutment	20°											
			45°											
	AstraTech EV	Fixture	3.0		1.27mm		—	P-DYNDRIVER						
			3.6											
			4.2											
			4.8											
	UniAbutment EV (DVU)	Abutment	Straight											
Multi base(DVM)	Abutment	Straight												
Ankylos (DY)	Fixture	C/X		1.0mm		—	N/A							
<b>Biohorizons</b>	Internal (BHB)	Fixture	3.5		1.25mm		✓	PADM14						
			4.5											
			5.7											
	Tissue level(BHT)	Fixture	3.5											
			4.5											
			5.7											
<b>Dentium</b>	Superline (DTS)	Fixture	—		1.28mm		✓	P-DYNDRIVER						
	Screw Abutment (DTU)	Abutment	4.5											
			5.5											
<b>Implant Direct</b>	active system	Fixture	3.0		1.25mm		✓	P-DYNDRIVER						
			3.4											
	Legacy system	Fixture	3.0											
			3.5											
			4.5											
	Raplant system	Fixture	5.7											
			3.5											
			4.3											
	Swissplant system	Fixture	5.0											
			6.0											
3.7														
4.8														
<b>Megagen</b>	Anyridge (MA)	Fixture	One Platform		1.20mm		—	N/A						
<b>MIS</b>	C1	Fixture	NP		1.27mm		✓	MT-ELM10						
			SP											
			WP											
	V3	Fixture	NP											
			SP											
Seven	Fixture	NP												
			SP											
			WP											
	Multi unit (MIU)	Abutment	4.8											
<b>Neodent</b>	CM	Fixture	One Platform		0.9mm		✓	105.148						
<b>Neoss</b>	—	Fixture	One		—		✓	51183 (22mm)						

## Common screw/screw driver systems of AS

Brand	System	Level	Platform diameter	Type of screw head (Straight ONLY)		Pictures of Straight Screw Hand Driver	ANGLED SCREW DRIVER:	Original Custom angled screw driver ref code:	Picture of ANGLED SCREW DRIVER					
Nobel Biocare	Replace Select	Fixture	NP		N/A		✓	300853 (28mm)						
			RP											
			WP											
	Active/Conical	Fixture	6.0											
			3.0											
			NP											
	Branemark	Fixture	RP											
			WP											
	MUA	Abutment	NP											
	On1	Abutment	RP											
WP														
Osstem	TS	Fixture	Mini		1.20MM		✓	ATD12L (20mm)						
			Standard											
Southern	External Hex (BO)	Fixture	3.0		1.22mm		—	N/A						
			3.25											
			3.75/4.0											
			5.0											
			6.0											
	Deep_Conical	Fixture	3.0											
			3.5											
			4.0											
	OCT-XT (ITT) Tapered	Fixture	4.8							—		✓	N/A	
			6.5											
	Tri-Nex	Fixture	3.5							—		✓	N/A	
			4.3											
			5.0											
	internal_hex	Fixture	3.75							1.27mm		—	N/A	
			4.2											
5.0														
Straumann	Tissue Level	Fixture	NNC		N/A		✓	O46.790 (27mm)						
			RN											
			WN											
	Bone Level	Fixture	SC											
			NC											
Screw-Retained Abutments (SBU)	Abutment	RC												
		3.5												
BLX (SX)	Fixture	4.6												
		RB												
		WB												
Zimmer	Tapered Screw-Vent	Fixture	3.5		1.25mm		✓	ASCDT30 (30mm)						
			4.5											
			5.7											